

WB-4MS, Stirred water bath



WB-4MS with base BP-1 (on the bottom)



Water bath-thermostat **WB-4MS** is designed for chemical, pharmaceutical, medical and biological laboratory research.

WB-4MS provides increased temperature stabilization (up to 0.1°C) due to built-in magnetic stirrer (regulated speed 250–1000 rpm).

Easy set up, high temperature maintenance accuracy, compact size and attractive modern design make this water bath widely used.

Specifications:

Temperature setting range	+25°C ... +100°C
Temperature control range	+5°C above ambient ... +100°C
Setting resolution	0.1°C
Stability	±0.1°C
Uniformity (at +37°C)	±0.1°C
Speed regulation range	250–1000 rpm
Digital timer with sound alarm	1 min–96 hrs
Temperature display	LCD, 2 x 16 signs
Digital setting of temperature, time and mixing speed	
Plastic lid with stainless steel interior included	
Quiet operation	
Tank capacity	4 litres
Working volume	235 x 135 x 110 mm
Overall dimensions	345 x 275 x 235 mm
Weight, not more	3.6 kg
Power supply	220/240 V, 50/60 Hz
Heating power	700 W

Catalogue number:

WB-4MS with base BP-1	BS-010406-AAA
------------------------------	---------------

